

# Airway Suction Trainer

M85



A representation of human respiratory organs, designed for temporary airway suction and catheter insertion practice. Effective tool with life-like tactile qualities for training in airway suction for choking patients.

## Features

1. Realistic anatomy includes the mouth, tongue, oral and nasal pharynx, epiglottis, trachea and esophagus.
2. Catheterization pathways include oral, nasal or through the tracheostomy opening.
3. Life-size organs facilitate learning proper catheter insertion length.
4. The side window allows the passage of the catheter to be observed.
5. Realistically simulated sputum is sticky to the touch.
6. The model is washable and can be easily dismantled for cleansing.

## Training Skills

- Temporary airway suction (intraoral and intranasal and endotracheal)



Nasal Suction



Oral Suction



Tracheal Suction

## Set Includes

- 1 head & neck model with airway anatomy
- 1 bottle of lubricant
- 2 catheters
- 1 bottle of simulated sputum
- 1 metal tray
- 1 syringe
- 1 storage case

## Manikin Size

W45 x D23 x H15 cm

## Replacement Parts

- 11229-020  
1 replacement skins with trachea
- 11229-010  
1 simulated sputum

- 11229-040  
2 sample catheters
- 11229-050  
1 lubricant